



VIVE

ALTERNATING SEAT CUSHION

CSH1084

HCPCS: E2610

The Vive Alternating Seat Cushion relieves pain and pressure on thighs and lower back, while providing maximum comfort when seated in a wheelchair or seat. Controlled by a quiet, rechargeable air pump, our pneumatic cushion provides a static setting and alternating airflow modes for redistributing weight across the cushion preventing pressure sores and relieving tension in the lower back.

SPECIFICATIONS

Product Dimensions:

- Cushion: 45 x 45 x 9cm (17.7" x 17.7" x 3.54")
- Pump Case: 27.5 x 13.5 x 8.5cm (10.8" x 5.3" x 3.3")

Product Material(s):

- Cushion - PU and Nylon/TPU material (waterproof)

Care and cleaning:

- Cushion: Spray with disinfectant then wipe clean with a damp cloth. Dry immediately
- Pump: Spot clean with a dry cloth (keep dry)

What's included:

- Vive Alternating Seat Cushion
- Rechargeable air pump
- Charging cable
- Replacement air cell
- Repair kit
- Fabric carry bag for pump

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



SIT COMFORTABLY FOR HOURS

The Vive Alternating Seat Cushion redistributes your weight and relieves painful pressure caused by prolonged seating. The pneumatic cushion also aids in the prevention of pressure sores and ulcers for optimal comfort when seated on a chair, wheelchair, recliner and outdoor seating.

CUSTOMIZABLE CUSHIONING

Our air seat cushion is equipped with six alternating cells providing alternating airflow to support the seat, thighs and lower back. The quiet air pump allows you to choose a static setting with a timer or three alternating modes for customizable comfort that safely supports up to 220 pounds.



PERFECT FOR WHEELCHAIRS

Sized to fit standard and bariatric wheelchairs and transport chairs, our alternating seat cushion is perfect for wheelchairs for adults, office chairs and more. The cushion also includes a protective waterproof cover that is easy to clean.



RECHARGEABLE AIR PUMP

Providing up to eight hours of cushioned seating, the Vive Alternating Seat Cushion's pressure pump is powered by a long-lasting rechargeable battery. The air pump is equipped with two hooks that allow it to be positioned on the back of the wheelchair or seat.



(239) 220-5772



sales@vivehealth.com

vive
health

SR22-12R1